



NU SLING™

Extension Strap

Directions to attach the extension strap for use with a cast or ace wrap



1. Slide the cuff strap all the way out of the rectangular D-ring. Lay the cuff out flat so the padding is facing up toward the ceiling. (The collar can be separated from the cuff or remain attached while adding the extension strap.)



2. Open the extension strap and feed it through the rectangular D-ring at the end of the wrist cuff, then secure to itself. Be careful to use the D-ring that you just removed from the cuff strap from and NOT the D-ring on the 3" adjustable loop strap.



3. Re-thread the original cuff strap through the new rectangular D-ring. You now have an additional 3" to use on your cuff.






NuSling - For orthopedic wrist, arm and shoulder injuries

NuSling goes beyond static immobilization and supports your dynamic recovery goals with its ergonomic design, one-handed operation and quick cuff release option in the MAGNETIC version.

Padded Collar



Padded Cuff

- A** Collar is made of wicking material and has an ergonomic shape to evenly disperse pressure across your neck and shoulders
- B** Collar can easily be slipped on/off overhead with one hand without using the front snap, but it is available as an option if needed
- C** 3" Velcro adjustable loop to change height of cuff/forearm up or down for positional comfort
- D** Adjustable Velcro cuff strap for securing the wrist or forearm
- E** Magnetic connection for quick release of collar and cuff to therapeutically use the arm without removing the sling (**Only** in NuSling Magnetic)



activearmsupports

NUSLING™

Comfort
Standard Cuff



Magnetic
Magnetic Cuff



Use & Cleaning: NuSling is intended for orthopedic shoulder, elbow or wrist injuries that require protection and modified use after an injury. Hand wash only with a mild soap and warm water. Wring out excess water and hang to dry.

Warning: NuSling Magnet contains two MAGNETS. If you have a pacemaker, defibrillator, neurotransmitter or other device implanted, please consult your doctor before using. Wear only as tolerated or recommended by your therapist.

Manufacturing/Contact Information: info@ActiveArmSupports.com; Active Arm Supports, P.O. Box 248, Hygiene, CO 80533; p: 303.519.2428