




NuSling - For orthopedic wrist, arm and shoulder injuries

NuSling goes beyond static immobilization and supports your dynamic recovery goals with its ergonomic design, one-handed operation and quick cuff release option in the MAGNETIC version.

Padded Collar



Padded Cuff

- A** Collar is made of wicking material and has an ergonomic shape to evenly disperse pressure across your neck and shoulders
- B** Collar can easily be slipped on/off overhead with one hand without using the front snap, but it is available as an option if needed
- C** 3" Velcro adjustable loop to change height of cuff/forearm up or down for positional comfort
- D** Adjustable Velcro cuff strap for securing the wrist or forearm
- E** Magnetic connection for quick release of collar and cuff to therapeutically use the arm without removing the sling (**Only** in NuSling Magnetic)



activearmsupports

NUSLING™

Comfort
Standard Cuff



Magnetic
Magnetic Cuff



Use & Cleaning: NuSling is intended for orthopedic shoulder, elbow or wrist injuries that require protection and modified use after an injury. Hand wash only with a mild soap and warm water. Wring out excess water and hang to dry.

Warning: NuSling Magnet contains two MAGNETS. If you have a pacemaker, defibrillator, neurotransmitter or other device implanted, please consult your doctor before using. Wear only as tolerated or recommended by your therapist.

Manufacturing/Contact Information: info@ActiveArmSupports.com; Active Arm Supports, P.O. Box 248, Hygiene, CO 80533; p: 303.519.2428



Directions for Use:



1. The NuSling will arrive fully assembled. Center the collar around your neck either by pulling it overhead or by using the front snap. Unhook the Velcro cuff strap and slide the involved hand/forearm through the open cuff. Move the cuff along your forearm to achieve a comfortable position, pull the strap down and resecure.

2. To move the cuff/forearm up or down, separate the Velcro pieces of the 3" adjustable loop and rotate to achieve the desired height. Press the Velcro adjustable loop together to hold the comfortable position. Small adjustments of the forearm up or down throughout the day can reduce aches and pain in the affected limb.

For NuSling Magnetic:

Easily separate the 2-pieces for functional use of the affected arm without removing the collar or cuff. Locate the small Velcro tab near the top of the of the 3" loop on the cuff. Use your UNaffected hand to pull the tab to unhook the small section of Velcro. Then hold the cuff/arm and push down and away from the collar to separate the magnets. Reverse the steps to reconnect the cuff and collar.