




LuxArm



Cleaning the Core Parts:

The belt pads and top cover of the forearm may be removed. The core belt and core forearm parts that remain are made of plastics and aluminum. These may be wiped down with an intermediate level disinfectant for sanitation purposes. Covers and pads may either be washed or replaced.

Directions for Replacement Pads/Cover:



Replacement Belt Pad Set
(1 Set = front pad, back pad and 2 hip pads)

The replacement belt pads come in a set of 4 pads that can be machine washed or replaced when worn out. Please connect any Velcro parts, wash in cold water, gentle cycle and tumble dry on low heat. The 2 hip pads slip over the corresponding hip templates on the core belt. The back pad is identified by the long Velcro strip. It is placed around the overlapping belt parts and button snaps in the back with the longer flap closing upward to connect with the shorter flap. The front pad is placed behind the strap with with ratchet buckle while the belt is open and the Velcro tab is closed around the belt to keep it from falling off. When it is time to close the belt the ratchet strap goes into the ratchet buckle while the plastic belt slides into the pocket of the front belt pad. Please be sure the padded aspect for each pad is on the inside where it will protect the skin and provide comfort.



Replacement Forearm Cover and Straps
(Thumb strap and extra finger strap not pictured)

The replacement forearm cover and straps can be machine washed or replaced when worn out. Please connect any Velcro parts, wash in cold water, gentle cycle and tumble dry on low heat. The cover slides on the easiest when the core forearm is contoured near the elbow and the wrist is bent into slight extension. Place the finger end of the forearm core into the cover then pull it over the contoured end. Be sure the end tab covers the core, but does not block the pin on the underside of the forearm core. The zipper tab on the cover has been removed to prevent it from interfering with the connection joint, so use a paperclip or small hook to pull the zipper closed on the underside of the forearm.